Stellest®

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H.A.L.T, Highly Aspherical Lenslet Target

* Two-year prospective, controlled, randomized, double-masked clinical trial results on 54 myopic children wearing Essilor® Stellest® lenses compared to 50 myopic children wearing single vision lenses in Wenzhou China.

Results based on 32 children from the Test Group wearing Essilor® Stellest® lenses at least 12 hours per day every day for two consecutive years

† Compared to combined single vision lens groups SVL(0M-24M)-SVL(24M-36M)

‡ Compared to the 60-month progression of the extrapolated control group (predicted average annual decrease in SER by 97%, Smotherman C, et al. IOVS 2023;64:ARVO E-Abstract 811 and predicted average annual decrease in AL by 15%.

§ Shamp W, et al. IOVS 2022;63:ARVO E-Abstract A011)

¶ One-year results from the two-year prospective, controlled, randomized, double-masked clinical trial results on 54 myopic children wearing Essilor® Stellest® lenses compared to 50 myopic children wearing single vision lenses.

‡ Subject to availability and launch timelines

By absorption. Additional UV back side reflection reduction when combined with Crizal® coating

1. Parssinen O, Kauppinen M, Viljanen A. The progression of myopia from its onset at age 8-12 to adulthood and the influence of heredity and external factors on myopic progression. A 23-year follow-up study. Acta Ophthalmol. 2014;92:730-739

2. Chua SY, Sabanayagam C, Cheung YB, Chia A, Valenzuela RK, Tan D, Wong TY, Cheng CY, Saw SM. Age of onset of myopia predicts risk of high myopia in later childhood in myopic Singapore children. Ophthalmic Physiol Opt. 2016;36:388-394

3. Matsumura S, Lanza C, Htoon HM, Brennan N, Tan C-S, Kathirani B, Chia A, Tan D, Sabanayagam C, Saw S-M. Annual Myopia Progression and Subsequent 2-Year Myopia Progression in Singaporean Children. Trans Vis Sci Tech. 2020;9:12-12

4. Donovan L, Sankaridurg P, Ho A, Naduvilath T, Smith EL, Holden BA. Myopia progression rates in urban children wearing single-vision spectacles. Optom Vis Sci. 2012;89:27-32

5. Jones LA, Sinnott LT, Multi DQ, Mitchell GL, Moeschberger ML, Zadnik K. Parental History of Myopia, Sports and Outdoor Activities, and Future Myopia. Invest Ophthalmol Vis Sci. 2007;48:3524-3532

6. Rose KA, Morgan IG, Ip J, Kiffey A, Huynh S, Smith W, Mitchell P. Outdoor Activity Reduces the Prevalence of Myopia in Children. Ophthalmol. 2008;115:1279-1285

7. Xiong S, Sankaridurg P, Naduvilath T, Zang J, Zou H, Zhu J, Lv M, He X, Xu X. Time spent in outdoor activities in relation to myopia prevention and control: a meta-analysis and systematic review. Acta Ophthalmol. 2017;95:551-566

8. Li SM, Li SY, Kang MT, Zhou Y, Liu LR, Li H, Wang YP, Zhan SY, Gopinath B, Mitchell P, Wang N, Anyang Childhood Eye Study G. Near Work Related Parameters and Myopia in Chinese Children: the Anyang Childhood Eye Study. PLoS One. 2015;10:e0134514

9. Zadnik K, Sinnott LT, Cotter SA, Jones-Jordan LA, Kleinstein RN, Manny RE, Twelker JD, Multi DQ. Collaborative Longitudinal Evaluation of E. Refractive Error Study G. Prediction of Juvenile-Onset Myopia. JAMA Ophthalmol. 2015;133:683-689

10. Wong YL, Li X, Huang Y, Yuan Y, Ye Y, Lim EW, Yang A, Spiegel DP, Drobe B, Bao J, Chen H. Eye growth pattern of myopic children wearing spectacle lenses with aspherical lenses compared with non-myopic children. Ophthalmic Physiol Opt. 2024;44(1):206-213

11. Bao J, Huang Y, Li X, Yang A, Zhou F, Wu J, Wang C, Li Y, Lim EW, Spiegel DP, Drobe B, Chen, H. Spectacle lenses with aspherical lenslets for myopia control vs single-vision spectacle lenses: a randomized clinical trial. JAMA ophthalmology. 2022;140(5):472-8

12. Li X, Huang Y, Yin Z, Liu C, Zhang S, Yang A, Drobe B, Chen H, Bao J. Myopia Control Efficacy of Spectacle Lenses with Aspherical Lenslets: Results of a 3-year Follow-up Study. Am J Ophthalmol. 2023;253:160-168

13. Drobe B, Xue L, Huang Y, Lim EW, Bao J. Spectacle lenses with highly aspherical lenslets for myopia control: 4-year clinical trial results. Investigative Ophthalmology & Visual Science. 2023;64(8):4162. Presented at the 2023 ARVO Annual Meeting, New Orleans

14. Li X, Huang Y, Liu C, Yin Z, Cui Z, Lim EW, Drobe B, Chen H, Bao J. Myopia control efficacy of spectacle lenses with highly aspherical lenslets: results of a 5-year follow-up study. Investigative Ophthalmology & Visual Science. 2024 Jun 17;65(7):31-17

15. Bao J, Yang A, Huang Y, Li X, Pan Y, Ding C, Lim EW, Zheng J, Spiegel DP, Drobe B, Lu F, Chen H. One-year myopia control efficacy of spectacle lenses with aspherical lenses. Br J Ophthalmol. 2022;106(8):1171-1176

©Essilor International - All Rights reserved - Do not copy or distribute - September 2024

This material is meant for the training of eye care professionals only, and is not to be reproduced or disclosed, whether with modifications or not, to any third parties or to the general public, without EssilorLuxottica's prior written consent.

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©Essilor® International - All Rights reserved - Do not copy or distribute - September 2024